

RTP4-receptor (chemosensory) transporter protein 4 Gene

Reference: PTXBC013161

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IFRG28; Z3CXXC4; receptor-transporting protein 4; 28 kDa interferon-responsive protein; 28kD interferon responsive protein; 3CxxC-type zinc finger protein 4; receptor (chemosensory) transporter protein 4; zinc finger, 3CxxC-type 4; receptor transporter pr		

For research use only.

GENE DESCRIPTION	receptor (chemosensory) transporter protein 4
GENE ID	64108
ORIGIN SPECIES	Human
ACCESSIONS	BC013161
VECTOR	pDONR223
NCBI SYMBOL	RTP4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtttagattctggaactgggagcagacattcaagaactcaagaggcaaaacccggggccacatggagcgtgaagttggatggcaacctcagctagactgcctggctcaagggtggaagcaatacaacagagagcattggctgggtccgtctcctgccagegaa
gttgggctccgccaagtgcagattctgtgccacacgtactgggagcactggacatccagggtcaggtgcgtatgaggctctttggccaaagggtccagaagtgcctcgtgtcccaatagagatgcctgagttctcctcggatagcaccatgaggattctgagcaacctgggtcagcat
atactgaagaatactatggaaatggcatgaggaagtctccagaatgccagtaatcctggaaggtgcctggaaggatcccatgacacagccaattgtgaggcatgcacttgggcataatgtgacagggccttaaaaagctacatgacaagccgtccaaatccctactccccacataaa
gactgggaattcctcacctggaattgtgtgtgtacctcgcaaccaagccaagaaccagtcagatgaggcaaaagggtgaagggtggtatgagaaattaggccagtcgagaccagatccactgaacatctgtgtctttatgtgtctgtattattgtagtcaaatgctttac
atcagaatga