

HLA-DPB1-major histocompatibility complex, class II, DP beta 1 Gene

Reference: PTXBC013184

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HLA-DP; HLA-DP1B; HLA-DPB; HLA class II histocompatibility antigen, DP beta 1 chain; HLA class II histocompatibility antigen, DP(W4) beta chain; HLA-DP histocompatibility type, beta-1 subunit; MHC HLA DPB1; MHC class II HLA-DP-beta-1; MHC class II antigen		

For research use only.

GENE DESCRIPTION	major histocompatibility complex, class II, DP beta 1
GENE ID	3115
ORIGIN SPECIES	Human
ACCESSIONS	BC013184
VECTOR	pDONR223
NCBI SYMBOL	HLA-DPB1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatgttctgcaggtttctgcggccccggacagtgggtctgacggcgttactgatggtgctgctcacatctgtgtccaggcaggccactccagagaattacgtgtaccaggacggcagggaatgtacgcgtttaatgggacacagcgttctctggagagatacatctacaacc
gggaggagtacgcgccttcgacagcgcgttggggaggttccggcggtgacggagctggggcgccctgctgcggagtagtgaacagccagaaggacatcctggaggagaagcgggcaatgcccggacagggatgcagacacaactacagctggacgaggccgtgacct
gcagcggcagtcacgctaaggtgaacgtttccctccaagaaggggccctgcagcaccacaacgtctgtctgccacgtgacagatttctaccaggcagcattcaagtcctgatggttcctgaatggacaggaggaaacagctggggtgctgtccaccaacgtgatccgaatgg
agactggacctccagatcctgtgtgctggaatgacccccagcaggagagcgtctacatctgccaaagtggagcacaccagcctggacagtcctgtcaccgtggagtggaggcacagtcgtattctgccagagtaagacattgacggagctggggccttctgctgggctca
tcactgtggagtgggcatcttcacagaggagcaagaagtcaacgaggtctgataa