

SLC25A38-solute carrier family 25, member 38 Gene

Reference: PTXBC013194

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SIDBA2; solute carrier family 25 member 38; appoptosin		
For research use only.			

GENE DESCRIPTION	solute carrier family 25, member 38
GENE ID	54977
ORIGIN SPECIES	Human
ACCESSIONS	BC013194
VECTOR	pDONR223
NCBI SYMBOL	SLC25A38
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgattcagaactcagtcctcgtcgtcgaacccaagatgtcggagacacgggtggaacgcttatgttacatcgggtgatcaaggcttctctgtgtggctccatcagtgaggacgtcctacccctctttccaacctctggatctccttaaacacgcctacaacccctcagccctcagat
catgggtctagacgtgttgggatgttgctgtactctgaaggtgttcgcacggagagctctttgggctttggaaggatgtccctccattgtgagatgttccctggcgttggaaatcactttggcactctctactctttgaagcagtatttctcgaggccatcccccaaccgcccctgga
gtcagtcagctcgtgggtgtgctcgtctgtgtcaggggtctgtatgtcacctatcactgtatcaagacgcgctatgagagtggaataatgctatgagagtatctacgtcgcctgaggagcatctacacagtgaggggcaccgggctctcagtgccctgacagcaactctct
tcgagatgcgccctctcaggaaatcactgatgtttacaaccagacaaaaatagtgccctatgaccaggtggatgaaccccttattctattacaaatttcagctgtgggatatttctgtgtattctggcctcactggtaactcaacctgcggatgttatcaaaactcatatgcagcttatcca
ctgaagttcaatggattggccaagcagtgacacttatttcaaaagactatggactacgtggcttctccaaggtggcatccccgagccctccgcagaactctaattggcagcaatggcgtggacggtgtatgaagagatgatggccaagatgggctgaagtcctga