

ATP5SL-ATP5S-like Gene

Reference: PTXBC013323

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ATP synthase subunit s-like protein; ATP5S like		
For research use only.			

GENE DESCRIPTION	ATP5S-like
GENE ID	55101
ORIGIN SPECIES	Human
ACCESSIONS	BC013323
VECTOR	pDONR223
NCBI SYMBOL	ATP5SL
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggctccctggcgctccctggcgctggcgcccccatgtggaatggcggtatcaggggcacatcgcctgggtgcgcagtggtggccagaggcgaatcagaagaagaaaggacaatactccagttcctgaccaactatttctacgatgtggaggctctgaggattactgtc
caaaggagatgtacaagggtcatgagaaaaatcggtttcagagacaaggatggatcagggccagataagtatggccatttctctcaggagttctggaatttctgtgaagtgcctgtcgaagctgtggatgccgtgactgtgacatcaactacgaggcctggataacctccgaacctcgc
aggctggacatctcggacctccctggcggtgtccaacctggcctcactcagatattgtggaggagatgtgcccaattgcgaggtgtgggagtcgactgggctgaggcctga