



Reference: PTXBC013325

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD365; HAVCR; HAVCR-1; KIM1; TIM; TIM-1; TIM1; TIMD-1; TIMD1; hepatitis A virus cellular receptor 1; T cell immunoglobulin domain and mucin domain protein 1; T-cell immunoglobulin mucin family member 1; T-cell immunoglobulin mucin receptor 1; T-cell membran		

GENE DESCRIPTION	hepatitis A virus cellular receptor 1
GENE ID	26762
ORIGIN SPECIES	Human
ACCESSIONS	BC013325
VECTOR	pDONR223
NCBI SYMBOL	HAVCR1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]