

GAS2-growth arrest-specific 2 Gene

Reference: PTXBC013326

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	growth arrest-specific protein 2; GAS-2; growth arrest specific 2		
For research use only.			

GENE DESCRIPTION	growth arrest-specific 2
GENE ID	2620
ORIGIN SPECIES	Human
ACCESSIONS	BC013326
VECTOR	pDONR223
NCBI SYMBOL	GAS2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgcactgctctgagcccaaggtacgagtgccctgctctgatgcacagtagcccaagacatgaagctaatgtctaccaatgaaagaagatctggccttggttaaccaatctattaggaaggagattacagcagaacttttatgagaagttgg
acaatggtgctctgctctgcaactgcagaaactatgcaggagaaattcaaggagagcatggatgctaacaagcccaagaagaatctaccgttgaagaagatcccatgcaaaaccagtgacccctcgggctcctttttgccagagacaatacagcaaatcttatctctggtgccgagattt
aggggtggatgaaacgtgtctatttgaatcggaaggttgggatttggaggccagtga