

TCAP-titin-cap (telethonin) Gene

Reference: PTXBC013330

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CMD1N; CMH25; LGMD2G; T-cap; TELE; 19 kDa sarcomeric protein; limb girdle muscular dystrophy 2G (autosomal recessive); teneurin C-terminal associated peptide; titin cap protein; titin-cap		

For research use only.

GENE DESCRIPTION	titin-cap (telethonin)
GENE ID	8557
ORIGIN SPECIES	Human
ACCESSIONS	BC013330
VECTOR	pDONR223
NCBI SYMBOL	TCAP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctacctcagagctgagctgcgaggtgtcggaggagaactgtgagcgcgggagggccttctgggcagaatggaaggatctgacactgtccacacggccgaggagggtgtccctgcatgaggaggacaccagagacatgagacctaccaccagcaggggcagtgccag
gtgctgggtgcagcgtcgccctggctgatgatcggtatgggcatccacggcgtgggctgcaggagtaccagctgcccaccacggggtactgcccgtgccatctcaccctgccaagatgggcgcaccaaggaggagcgtgaggacaccccatcagcttcaggagctgctg
gcgctggagacagccctgggtggcaggtgtgtggaccgccaggaggtggctgagatcacaagcagctgcccctgtgtgctgtcagcaagccgggtgcccttcgtcgtccctgtcccgtccatgtcccaggaagcacagagaggctga