

RAB11A-RAB11A, member RAS oncogene family Gene

Reference: PTXBC013348

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RAB11A, member RAS oncogene family; YL8; ras-related protein Rab-11A; RAB11A, member oncogene family; rab-11		

For research use only.

GENE DESCRIPTION	RAB11A, member RAS oncogene family
GENE ID	8766
ORIGIN SPECIES	Human
ACCESSIONS	BC013348
VECTOR	pDONR223
NCBI SYMBOL	RAB11A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggcacccgcgacgacgagtagcactacctctttaagttgtccttattggagattctggtgttggaagagtaatctctctgatttactgaaatgagtttaactctggaagcaagagcaccattggagtagagtttgcaacaagaagcattccaggtgatggaaaaacaataaaggc
acagatatgggacacagcagggcaagagcgatcgaactataacatcagcatattatctgtggagctgtagtgccctattgtttatgacattgctaacaatcctacatatgaaaatgtagagcgatggctgaaagaactgagagatcctgctgataacattgttatcatgcttgggcaa
taagagtgatctacgtcatctcagggcagttcctacagatgaagcaagagcctttgcagaaaaagaatggtttgtcattcattgaaacttcggccctagactctacaatgtagaagctgcttttcagacaattttaacagagattaccgcattgtttcagaagcaaatgtcagacagacgcgaa
aatgacatgtctccaagcaaatgtggttctattcatgttccaccaaccactgaaaacaagcgaaggtgcagtgctgtcagaacatctaa