



Reference: PTXBC013352

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HTF9C; tRNA (uracil-5-)-methyltransferase homolog A; TRM2 tRNA methyltransferase 2 homolog A; hpaII tiny fragments locus 9c protein; tRNA methyltransferase 2 homolog A		

For research use only.

GENE DESCRIPTION	TRM2 tRNA methyltransferase 2 homolog A (S. cerevisiae)
GENE ID	27037
ORIGIN SPECIES	Human
ACCESSIONS	BC013352
VECTOR	pDONR223
NCBI SYMBOL	TRMT2A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]