



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

FOXS1-forkhead box S1 Gene

Reference: PTXBC013408

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FKHL18; FREAC10; forkhead box protein S1; FREAC-10; forkhead-like 18 protein; forkhead-related activator 10; forkhead-related transcription factor 10; forkhead-related transcription factor FREAC-10; forkhead box S1		

For research use only.

GENE DESCRIPTION	forkhead box S1
GENE ID	2307
ORIGIN SPECIES	Human
ACCESSIONS	BC013408
VECTOR	pDONR223
NCBI SYMBOL	FOXS1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr