

TMEM208-transmembrane protein 208 Gene

Reference: PTXBC013412

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HSPC171; transmembrane protein 208		
For research use only.			

GENE DESCRIPTION	transmembrane protein 208
GENE ID	29100
ORIGIN SPECIES	Human
ACCESSIONS	BC013412
VECTOR	pDONR223
NCBI SYMBOL	TMEM208
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgcccaagggcaagtgggcacgagaggggaagaagcagatattgaagagaacagagagactctgaagtctacctgcggatcatactggggccaatgccattactgccttgtagcgttggtctctttactatctgcctcatttgggcctggtggccctgggctttagtctggc
agtgtatggggccagctaccactctatgagctcgaatggcacgagcagcgttctctgagatgagggccctgatggatggatggacccaacatggagcagggcatggcagagcaccttaaggatgtatctactgacagccatgtgcagggtgctcagctgcttctctctatgtctg
gtccttctggctctgagcagggccggttacctctgtgggtgaatgtctgggcccctggttctactgcagacagtgccacccagcaccagagcacaatgagaacggcagcgccgacaggagcggcgagatgaagcgttatag