

SELM-selenoprotein M Gene

Reference: PTXBC013421

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	selenoprotein SelM; SELM; SEPM; selenoprotein M		

For research use only.

GENE DESCRIPTION	selenoprotein M
GENE ID	140606
ORIGIN SPECIES	Human
ACCESSIONS	BC013421
VECTOR	pDONR223
NCBI SYMBOL	SELM
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagcctcgtgtgccctccgctggcgtgctgctgtctcgcggcgctgtggccccagccacagccgcaactgcctaccggcggaactggaaccgtctgagcggcctaaccgcgcccggttagagacctgcgggggatgacagctgaaccgcctaaggaggtgaagcttcg
tcacgcaggacattccattatcacaacctgggtgatgaacacctccctggggccgacctgagctgctgctgctggcgccgctacgaggaactagagcgcacccactcagtgaaatgacccgcgaagagatcaatgcgctagtgcaggagctcgcttctaccgcaaggcgccg
cccgcgcgcagggtgcccccgagtagctgtgtggcgcccggaagccccagaggaaactcggaccacgctgacctgtag