

SPIN1-spindlin 1 Gene

Reference: PTXBC013571

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SPIN; TDRD24; spindlin-1; ovarian cancer-related protein; spindlin1; spindlin 1		
For research use only.			

GENE DESCRIPTION	spindlin 1
GENE ID	10927
ORIGIN SPECIES	Human
ACCESSIONS	BC013571
VECTOR	pDONR223
NCBI SYMBOL	SPIN1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagacccattcggaagacacctggccagcggtccagagctgatgcaggccatgctggagtatctgccaacatgatgaagaaggagacatccacaaaaacatcggagcagtggtgggtccgagcaaacctgtttccagccccggcggaacatcgtaggctgcaggattcag
catgggtggaagaggggaatggccctgttacccagtggaagggaacgcttctggaccaggtgctgtaaatcctctttgtatctataaaatagcagtgattgactgtgtttatggactagaacttaataaagatgaaagagtttctgcgctgaagtcctccctgatagattgcgacatc
gaatcagcgatgcacactggcagacacaatgattggcaagcagtggaacatattgttgagacagagagtggttctaaagatgagtgagggggaatggtcttagcacgtgcacgtgcatgaacacatgggtttacattacatgagaaagacccctgtctgtacatgtaccaactcttagat
gattacaagaaggcgacctgcattatgcctgattccaatgattcacctccagcagaagggaaccaggagaagttgtggacagcctggtaggcaacaagtggaatatgcacaagaagatgggtcgaaggactggcatggtcatcatcaagtagaagccaagccctccgtctatt
tcatcaagtttgatgatgattccatatttatgtctacgatttggtgaaacatcctag