

CLDN11-claudin 11 Gene

Reference: PTXBC013577

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	OSP; OTM; claudin-11; oligodendrocyte transmembrane protein; oligodendrocyte-specific protein; claudin 11		

For research use only.

GENE DESCRIPTION	claudin 11
GENE ID	5010
ORIGIN SPECIES	Human
ACCESSIONS	BC013577
VECTOR	pDONR223
NCBI SYMBOL	CLDN11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtggccacgtgcctgcagggtgggtcttcacgagcttcgtgggtggatcgggtcctcgtgaccacctccaccaatgactgggtggtgacgtcggtacacctcccacgtgccgaagctggatgagctgggtccaaggggtgtgggtggcgtgcgtcatggccac
gggggtgtaccactgcaagccctgtgtgacatcctcatctgcgggtacgtgcaggcctgcgcgcctgatgattgctgcctcgggtcgtggtcgtgcccgtatcttactgctgactgtttcctcctgcatccggatgggcccaggagcccggtgtggctaagtacaggcgggcc
cagctggctgtgtgtgtcattctgctgctctgcgccttgttgcacacatctggttcctgtgtgcgcccaccgtgagaccacatctgagcttggctactccctgtatgcaggctggatggtgctgtgtgcctcgtgggtggctgtgtcatcctctgctgcgtggagatgcc
aggccttgggtgaaaccgtttctactacactgcgggtctagctccccgactcatgcgaagagtggccacgtataa