



Reference: PTXBC013586

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AGE-R2; G19P1; GIIB; PCLD; PCLD1; PKCSH; PLD1; VASAP-60; glucosidase 2 subunit beta; AGE-binding receptor 2; advanced glycation end-product receptor 2; glucosidase II subunit beta; hepatocystin; protein kinase C substrate 60.1 kDa protein heavy chain; pro		

GENE DESCRIPTION	protein kinase C substrate 80K-H
GENE ID	5589
ORIGIN SPECIES	Human
ACCESSIONS	BC013586
VECTOR	pDONR223
NCBI SYMBOL	PRKCSH
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]