

C10orf35-chromosome 10 open reading frame 35 Gene

Reference: PTXBC013587

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C10orf35; chromosome 10 open reading frame 35		
For research use only.			

GENE DESCRIPTION	chromosome 10 open reading frame 35
GENE ID	219738
ORIGIN SPECIES	Human
ACCESSIONS	BC013587
VECTOR	pDONR223
NCBI SYMBOL	C10orf35
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgcgcatcttgccaatggggaatcgtgcaggatgacgaccccgagtgaggaccactaccagccaccaagagtagcattcctcgacagagcttctcaataggggcatgggtgctccccaggggtcctggccccgccagcagcaggcagggtgccaggctgggtgct
gctcagtcctcccttcaatgacctcaaccggcagctgggaacatgggcttccgagtgagcatctcggaacacatgctgtggagccggtgacctccctctcctcttctgctcatgatgcttggtgttcgtggcctcctcctgggtggccttgtctacctggtgtccacctgagtcagcgg
tga