

UNC84A-unc-84 homolog A (C. elegans) Gene

Reference: PTXBC013613

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UNC84A; SUN domain-containing protein 1; Sad1 unc-84 domain protein 1; protein unc-84 homolog A; sad1/unc-84 protein-like 1; unc-84 homolog A; Sad1 and UNC84 domain containing 1		

For research use only.

GENE DESCRIPTION	unc-84 homolog A (C. elegans)
GENE ID	23353
ORIGIN SPECIES	Human
ACCESSIONS	BC013613
VECTOR	pDONR223
NCBI SYMBOL	UNC84A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatttttctcggcttcacatgtacagtcctcccagtggtgcccggagaacacgggctacacgtatgcgtcagttccagctattcttcagatgctctggattttgagacggagacacaaattggacctgtattgattctccacggatgtcccgcgtagtttgcgcctggccacgacagcat
gcacctgggggatggtgagcgtgtgggtgcccacagcggcaccagcagcgtgtctccctgaagaaccgagcggccagaaacaaacagcgcagaaacacaaatcagcttttagtatcaaccacgtgtcaaggcaggtcacgtctcttggcgtcagccacggcgccact
gtcagcctgcaggatgctgtgactcgacggcctctgtattggacgagtcgttgattcgtgaacagacacagtgaccactctctgggtcttgatgatggtgatcttaagggtgaaataaaagctgccattcagggaaacggggatgtggagccgccgccaccgcgcacaac
ggcttctctgcagcaactgcagcatgctgtccgagcgcaaggacgtgctcagcggcgaccccgcccccggcccggtgtcagagatttatctagggacaggaaatcaaaatgtaagtctcagtcctttaaactcagaaaaagggtgttttccaatttaatttcttctgtaagtct
cagtgctgcactattgtcttgagacttaaaattatcccttga