



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

JTV1-JTV1 gene Gene

Reference: PTXBC013630

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

JTV1; JTV-1; P38; aminoacyl tRNA synthase complex-interacting multifunctional protein 2; ARS-interacting multi-functional protein 2; multisynthase complex auxiliary component p38; multisynthetase complex auxiliary component p38; protein JTV-1; aminoacyl t

For research use only.

GENE DESCRIPTION	JTV1 gene
GENE ID	7965
ORIGIN SPECIES	Human
ACCESSIONS	BC013630
VECTOR	pDONR223
NCBI SYMBOL	JTV1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccgatgtaccaggctaaagccctatcacggggggggcgctctccgtgtggagctctccacctgcatgtaccggctcccaacgtgcacggcaggagctacggccacgcgcggggcgctggccacgtgcaggaaagactaacctgtctctgcaagctcttgatcccgcaag
atgataattaaacgctgtcatgtgttgaaagacagctgattgtggctctccaagatgtatcaaacacacagatgcacgtgtgattgataaacacataatcaacggcgatgagccacagacttaacacacaaatgcgcgtgacttgatgaattcagctgggaagattacggggcgctctgaaga
catcgtgttaacgcgaacacccctggctccctccctctccctctctctgctgacacgctctctctgtgacatgcctcagctggctctgccacgtgcacacgacacctctcctcagcagcagctgcctgaaacctctcgaagctgcttgcttgaggagacagataaaaaaacacggccgcgaagactatc
agctcggattcacttaatttgggaagatgtccgaagacgcagatgaaatcagacaccagacagatgtgccccatcgaaggcgaagggaacattgcacgtttctgttctctctgttggccagaagcataatgctgtcaacgcaaccttatagatagctgggtagatattcgattttcagtt
aaaagaggaagacagtaagaaaaagggcctgttttccgctccatgaactctgctcttgggaagacgccttgcctctgggaatgaactaccgtgacacagctgctgtgtctgtactccacagatcggagcgtcagctgtgacatgtcccgacccaatgtgcagacgtggatgag
gtcttctgaaacctgcctctttaacacgcgcctccaacgtcttaagtgta

More info: Please contact us for any question: contact@proteogenix.fr