

TXN2-thioredoxin 2 Gene

Reference: PTXBC013726

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	COXPD29; MT-TRX; MTRX; TRX2; thioredoxin, mitochondrial; thioredoxin 2		
For research use only.			

GENE DESCRIPTION	thioredoxin 2
GENE ID	25828
ORIGIN SPECIES	Human
ACCESSIONS	BC013726
VECTOR	pDONR223
NCBI SYMBOL	TXN2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctcagcgaactctctgaggaggtctggtcctctcatctccagggaagccctcagggtcagtggtccaccctcacttccagagccctgcagacccacacatgcagtcctggtggtcgtgactgtaacaccaaccagcccgacaatatacaccacgaggatctctgacaac
ctttaataccagatggacgtgactttcaagaccgagtggtcaacagtgagacaccagtggtgtggtgattccacgcacagtggtgtgacccgtgcaagatcctggggccgaggttagagaagatggtgccaagcagcacgggaaggtggtgatgccaaggtggatattgatgacca
cacagacctgccaatgagtgaggtgtcagcgggtgcccactgtgctggtccatgaagaatggggacgtggtggacaagttgtgggcatcaaggatgaggatcagttggagccctcctgaagaagctgattggctga