

ZNF501-zinc finger protein 501 Gene

Reference: PTXBC013762

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ZNF; ZNF52; zinc finger protein 501; zinc finger protein 52		
For research use only.			

GENE DESCRIPTION	zinc finger protein 501
GENE ID	115560
ORIGIN SPECIES	Human
ACCESSIONS	BC013762
VECTOR	pDONR223
NCBI SYMBOL	ZNF501
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaacatgggagagtttaacatgcagaagaacctcaaaagtgtatgaatgtgggaagttcttactcagagatcatcttaccagcaccagaggattcacagaggagagaagccctatgtgtcagtgatgtggaagttgttccgtaaacagtcacaaatctactcaacatctgagaat
tcataccggagagaaaccttataatgtaatgaatgtgagaagccttcaacaaaagcaattctgttcagcatctgagaattcactcaggagagaacctataatgcaatgaatgtggaaagcctttgtcagagcccatccctattaaacaccagcgaattcactcaggagaaacc
atataaattacagaatgtggcaaaccttcagtcagagcatatgccttactcgtcatcagagaagtcattctggagataaaccttttaagtgtatgaatgtggaaagccttttaacagagtcagtgctcagcatcagagaattcattcaggagagaagccctacacatgcactgaatg
tggtaaaccttcactcagaactctccctgttgaacatgaaaggactcacatggagagaacctttataagtgtagtgagtggtgaaaaaactttccgcaacaagcacacccttagtgagcattacagaattcactcaggagaaaaccttatgagtggttgatgtgggaatccttaggca
cagttcagcacttcttcacatcagagccttcactcgtggagagtaa