

C4orf42-chromosome 4 open reading frame 42 Gene

Reference: PTXBC013775

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C4orf42; CTBP1-AS1; CTBP1 antisense RNA 2 (head to head)		
For research use only.			

GENE DESCRIPTION	chromosome 4 open reading frame 42
GENE ID	92070
ORIGIN SPECIES	Human
ACCESSIONS	BC013775
VECTOR	pDONR223
NCBI SYMBOL	C4orf42
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaccaccgcccagggcagcgtttggcgtgtgtgtgcctgccaccgccctcctgtgcctggcctgccctgtgcctgcacacaggagggcctcatggcgtttgccctacacggatggcgtctcctggagctactggacagtcaccttgggggaatgccagagcatgggc
attaggtcccccgccagcctccgttgcacatgggctattttgtccatcgcgtgggacaacctagtattgggggaaaactcagtcactctaa