



Reference: PTXBC013781

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	104p; GCP3; Grip104; SPBC98; Spc98; Spc98p; gamma-tubulin complex component 3; gamma-ring complex protein 104 kDa; spindle pole body protein Spc98 homolog; tubulin gamma complex associated protein 3		

For research use only.

GENE DESCRIPTION	tubulin, gamma complex associated protein 3
GENE ID	10426
ORIGIN SPECIES	Human
ACCESSIONS	BC013781
VECTOR	pDONR223
NCBI SYMBOL	TUBGCP3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]