

TMEM185A-transmembrane protein 185A Gene

Reference: PTXBC013793

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CXorf13; FRAXF; ee3; transmembrane protein 185A; family with sequence similarity 11, member A; fragile site, folic acid type, rare, fra(X)(q28) F		

For research use only.

GENE DESCRIPTION	transmembrane protein 185A
GENE ID	84548
ORIGIN SPECIES	Human
ACCESSIONS	BC013793
VECTOR	pDONR223
NCBI SYMBOL	TMEM185A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccttctgtgcctgggtgctctactacattgtgtggcgtctgtttcttgcgtctatggatgtgaltgcggagcgcaggacacataacatggccctgagctggatgaccatcgtcgtgccctcttacatttgagattctgtggttcacaaactggatggccacaacgccttc
tcctgcatcccgatctttgtccccccttggctctcgttgatcacgctgatggcaaccacatttggacagaaggaggaaaccactggtgttggatccgcaaaagtctcagtttctgtgaaatctccatttctacgagaatatggaacatttctatgatccatcacgaagataatg
aagaaaccgaagagaccagttccggagcccctaaatcgcacccatgtttcgaagaaggccagggtgtcattaccagagccctgggaagtgtgtctccacctcccaattaaatatgaaatgccagattag