



Reference: PTXBC013885

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	4.1N; MRD11; band 4.1-like protein 1; neuron-type nonerythroid protein 4.1; neuronal protein 4.1; erythrocyte membrane protein band 4.1 like 1		

For research use only.

GENE DESCRIPTION	erythrocyte membrane protein band 4.1-like 1
GENE ID	2036
ORIGIN SPECIES	Human
ACCESSIONS	BC013885
VECTOR	pDONR223
NCBI SYMBOL	EPB41L1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]