

TTC26-tetratricopeptide repeat domain 26 Gene

Reference: PTXBC013912

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DYF13; IFT56; dyf-13; intraflagellar transport protein 56; TPR repeat protein 26; tetratricopeptide repeat protein 26; tetratricopeptide repeat domain 26		

For research use only.

GENE DESCRIPTION	tetratricopeptide repeat domain 26
GENE ID	79989
ORIGIN SPECIES	Human
ACCESSIONS	BC013912
VECTOR	pDONR223
NCBI SYMBOL	TTC26
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatgcttcaaggccaacacgtctgtaggcagagcgctacagcacactgacaaaagaagaagaagtaggaagattccaaactagaggagctcttcaaaaagagatttcaactgagctattaccctgttgagttcaaacgtcatgttggggaagaagaaggataactaattgttgattggatattgtgctttcacctgggtgactacaagagagctctggaggaatacgaataatgctcaaaaagaggaaaattgtaattctgaagctctgggtgaacctagcttgacactactcttcttgggatgtataaacaagctgaagcagctggatttaagctcaaaaagccgactcca
aaaccgcctcctctccacttggctcacaagtttaattgatgagaaaaaattgatgagctttcatcaaaatcttcaggatgtcacagaagatcaactcagtttggcctcaatccactatatgcgactcactaccaagaagctatagatatataagcgaatactgctagataacagggaatacttg
ccctaatgtttatgtggcctctgctactacaagttgattactatgatgtgtcgaagaagtttggcgtttaccttcagcaaatcctgatagaccatgcactcaatcttaagcctgaaccatttgcctttacaatggcagagcagctgaggcagaactcaaaagcttgatgacaatgct
tcttcactcttgaatttgcataaagaactcatcaggcacaatctggtaataacgatttccgccttga