



Reference: PTXBC013983

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HSAJ1454; TES-3; TICN3; testican-3; sparc/osteonectin, cwcv and kazal-like domains proteoglycan (testican) 3; SPARC/osteonectin, cwcv and kazal like domains proteoglycan 3		

For research use only.

GENE DESCRIPTION	sparc/osteonectin, cwcv and kazal-like domains proteoglycan (testican) 3
GENE ID	50859
ORIGIN SPECIES	Human
ACCESSIONS	BC013983
VECTOR	pDONR223
NCBI SYMBOL	SPOCK3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]