



Reference: PTXBC014021

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BT2.2; BTF2; BTN2.2; butyrophilin subfamily 2 member A2; butyrophilin 2		
For research use only.			

GENE DESCRIPTION	butyrophilin, subfamily 2, member A2
GENE ID	10385
ORIGIN SPECIES	Human
ACCESSIONS	BC014021
VECTOR	pDONR223
NCBI SYMBOL	BTN2A2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]