



Reference: PTXBC014026

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HHR23A; HR23A; UV excision repair protein RAD23 homolog A; RAD23, yeast homolog, A; RAD23 homolog A, nucleotide excision repair protein		

For research use only.

GENE DESCRIPTION	RAD23 homolog A (<i>S. cerevisiae</i>)
GENE ID	5886
ORIGIN SPECIES	Human
ACCESSIONS	BC014026
VECTOR	pDONR223
NCBI SYMBOL	RAD23A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]