



Reference: PTXBC014058

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	mitochondrial intramembrane-cleaving protease PARL; PRO2207; PSARL; PSARL1; PSENIP2; RHBDS1; presenilins-associated rhomboid-like protein, mitochondrial; rhomboid 7 homolog 1; presenilin associated rhomboid like		

For research use only.

GENE DESCRIPTION	presenilin associated, rhomboid-like
GENE ID	55486
ORIGIN SPECIES	Human
ACCESSIONS	BC014058
VECTOR	pDONR223
NCBI SYMBOL	PARL
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr