

RAP1A-RAP1A, member of RAS oncogene family Gene

Reference: PTXBC014086

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RAP1A, member of RAS oncogene family; C21KG; G-22K; KREV-1; KREV1; SMGP21; ras-related protein Rap-1A; GTP-binding protein smg p21A; Ras-related protein Krev-1		

For research use only.

GENE DESCRIPTION	RAP1A, member of RAS oncogene family
GENE ID	5906
ORIGIN SPECIES	Human
ACCESSIONS	BC014086
VECTOR	pDONR223
NCBI SYMBOL	RAP1A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcgtgagtacaagctagtggtccttgggtcaggagcgttgggaagtctgctctgacagttcagttgttcagggaattttgttgaataatgacccaacgatagaagattcctacagaaagcaagttgaagtcgattgccacagtgatgctcgaatcctggatactgcagggacaga
gcaatttacagcaatgagggatttggatatagaagaacggccaaggttttgactagtagtattctattacagctcagtcacgtttaacgacttacaggacgtgagggaacagattttacgggttaaggacacggaagatgttccaatgattttgttggcaataaattgacctggaagatgagcg
agtagttggcaagagcagggccagaatttagcaagacagtggtgtaactgtgctcttttagaattctctgcaagtcgaagatcaatgttaatgagatattttatgacctggtcagacagataataggaaaaccagtggaagaaagaagcctaaaaaagaatcatgtctgctctag