



Reference: PTXBC014100

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GCC1P; GCC88; GRIP and coiled-coil domain-containing protein 1; golgi coiled-coil 1; golgi coiled-coil protein 1; peripheral membrane Golgi protein; GRIP and coiled-coil domain containing 1		

For research use only.

GENE DESCRIPTION	GRIP and coiled-coil domain containing 1
GENE ID	79571
ORIGIN SPECIES	Human
ACCESSIONS	BC014100
VECTOR	pDONR223
NCBI SYMBOL	GCC1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]