



Reference: PTXBC014109

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	leucine-rich repeat-containing protein 45; leucine rich repeat containing 45		
For research use only.			

GENE DESCRIPTION	leucine rich repeat containing 45
GENE ID	201255
ORIGIN SPECIES	Human
ACCESSIONS	BC014109
VECTOR	pDONR223
NCBI SYMBOL	LRRC45
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr