



Reference: PTXBC014137

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein 7-60; zeta-type opioid receptor		
For research use only.			

GENE DESCRIPTION	opioid growth factor receptor
GENE ID	11054
ORIGIN SPECIES	Human
ACCESSIONS	BC014137
VECTOR	pDONR223
NCBI SYMBOL	OGFR
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr