



Reference: PTXBC014152

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CK-6C; CK-6E; CK6A; CK6C; CK6D; K6A; K6C; K6D; KRT6C; KRT6D; PC3; keratin, type II cytoskeletal 6C; cytokeratin 6A; cytokeratin 6C; cytokeratin 6D; keratin 6A, , type II; keratin 6A, type II; keratin, epidermal type II, K6A; keratin, type II cytoskeletal		

GENE DESCRIPTION	keratin 6A
GENE ID	3853
ORIGIN SPECIES	Human
ACCESSIONS	BC014152
VECTOR	pDONR223
NCBI SYMBOL	KRT6A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]