

BAX-BCL2-associated X protein Gene

Reference: PTXBC014175

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	apoptosis regulator BAX; BCL2L4; BCL2 associated X protein; BCL2-associated X protein omega; Baxdelta2G9; Baxdelta2G9omega; Baxdelta2omega; bcl-2-like protein 4; bcl2-L-4; BCL2 associated X, apoptosis regulator		

For research use only.

GENE DESCRIPTION	BCL2-associated X protein
GENE ID	581
ORIGIN SPECIES	Human
ACCESSIONS	BC014175
VECTOR	pDONR223
NCBI SYMBOL	BAX
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacgggtccggggagcagcccagaggcggggggccaccagctctgagcagatcatgaagacagggcccttttctcagggtttcatccagatcagcagggcggaatggggggggaggcaccgagctggccctggaccgggtgcctcaggatgcgtccaccaagaag
ctgagcgagtgctcaagcgcacgcggggacgaactggacagtaacatggagctgcagaggatgattgcccgtggacacagactcccccgagaggtcttttccgagtggcagctgacatgttttctgacggcaactcaactggggccgggttgcgcccttttctacttgcagcaa
actggtgctcaaggccctgtgcaccaaggtgcccgaactgacagaaccatcatggcctggacattggacttctccgggagcggctgttgggctggatccaagaccagggtgttgggacggcctcctctctacttgggacggccacgtggcagaccgtgaccatcttggcggg
agtgtccaccgcctcactcaccatctggaagaagatggcctga