

RAB30-RAB30, member RAS oncogene family Gene

Reference: PTXBC014213

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RAB30, member RAS oncogene family; ras-related protein Rab-30		
For research use only.			

GENE DESCRIPTION	RAB30, member RAS oncogene family
GENE ID	27314
ORIGIN SPECIES	Human
ACCESSIONS	BC014213
VECTOR	pDONR223
NCBI SYMBOL	RAB30
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagtatggaagattatgatttctgttcaaaattgtttaattggcaacgtgtgtgtgggaagacgtgcctcgtccgaagattcactcagggtcttttccccaggtcaaggagccacaattggagttgattttatgattaagacagtgagattaatgggtgaaaaagtaagctacagatctg
ggacacagcaggtcaagagagatttcggtccattaccagagttactaccgaagcgccaatgccttgatcctcacctatgacattacgtgtgaggaatccttccgttgccttcctgagtgctgaggagatagaacaatatgccagcaacaaggtcatcactgtgttagtgggcaacaagatt
gacctggctgaaaggagagaggtttccagcagcagctgaagaattctcagaagctcaggacatgattatctggagacctcagccaaggaatctgataatgtggagaaactcttcttgacttagcatgccgactcatcagtgaagccagacagaacacacttgaacaatgtatcctca
cccttacctggagaagggaagcatcagctatttgactgtgttaattcaactaa