



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

FAM113A-family with sequence similarity 113, member A Gene

Reference: PTXBC014247

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FAM113A; C20orf81; bA12M19.1; PC-esterase domain-containing protein 1A; family with sequence similarity 113, member A; sarcoma antigen NY-SAR-23; PC-esterase domain containing 1A		

For research use only.

GENE DESCRIPTION	family with sequence similarity 113, member A
GENE ID	64773
ORIGIN SPECIES	Human
ACCESSIONS	BC014247
VECTOR	pDONR223
NCBI SYMBOL	FAM113A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr