

PRR13-proline rich 13 Gene

Reference: PTXBC014257

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TXR1; proline-rich protein 13; taxane-resistance protein; proline rich 13		
For research use only.			

GENE DESCRIPTION	proline rich 13
GENE ID	54458
ORIGIN SPECIES	Human
ACCESSIONS	BC014257
VECTOR	pDONR223
NCBI SYMBOL	PRR13
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtggaatcccaatgccggcagccaggccaaatccatcccccaatattgggtgccctggaggttccaatcctgcccaaccacacattataatccaccttccccaggccccctgtcctctccccaggagctccccatggcaatccagcttccccagggtggccccctca
tcctgtgccacagccagggtatccagatgccaaaccgttgggtccctaccctcctccataccacccgctgccctggaaatccctcctgtgaatcccttggctcctggcatgttggaccagcagtgatagtagacaagaagatgcagaagaaaatgaagaaagctataaaagatgcac
aagcaccaaaagcaccacaagfaccacaagcatggcaagcattcctctctctcctcctctccagcagtgattctgactga