

RPL37A-ribosomal protein L37a Gene

Reference: PTXBC014262

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	L37A; 60S ribosomal protein L37a		

For research use only.

GENE DESCRIPTION	ribosomal protein L37a
GENE ID	6168
ORIGIN SPECIES	Human
ACCESSIONS	BC014262
VECTOR	pDONR223
NCBI SYMBOL	RPL37A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccaaacgtaccaagaaagtcgggatcgctgtaaacgggacccgctatggggcctccctccggaaaatggggaagaaatgaaatcagccagcagccaagtacacttgcctttctgtggcaaaccaagatgaagagacagactgtggggatctggcactgtggttcctgc
atgaagacagtgctggcgtgctggcgtgacgtacaataccacttccgctgacggtgtaagtccgcatcagaagactgaaggagtgaaagaccagtag