



Reference: PTXBC014271

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	END; HHT1; ORW1; endoglin; CD105 antigen		
For research use only.			

GENE DESCRIPTION	endoglin
GENE ID	2022
ORIGIN SPECIES	Human
ACCESSIONS	BC014271
VECTOR	pDONR223
NCBI SYMBOL	ENG
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]