

TIMP3-TIMP metalloproteinase inhibitor 3 Gene

Reference: PTXBC014277

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HSMRK222; K222; K222TA2; SFD; metalloproteinase inhibitor 3; MIG-5 protein; protein MIG-5; tissue inhibitor of metalloproteinases 3; TIMP metalloproteinase inhibitor 3		

For research use only.

GENE DESCRIPTION	TIMP metalloproteinase inhibitor 3
GENE ID	7078
ORIGIN SPECIES	Human
ACCESSIONS	BC014277
VECTOR	pDONR223
NCBI SYMBOL	TIMP3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaccccttggtcgggtcctcgtctcctgggcagctggagcctgggggactgggggcccggcgctgcacatgctgccagcccccaggagcctctgcaactccgacatcgtgatccggccaaggtggtgggaagaagctgtaaggaggggcccttcggcac
gctggtctacacatcaagcagatgaagatgaccgaggttcaccaagatgccccatgtgcagttacatccatcaggaagcttcgagagttctgtggtccttaagctggaggtcaacaagtaccagtacctgctgacaggctcgtctatgatggcaagatgtacacgggctgtgcaact
tcgtggagaggtgggaccagctcaccctctccagcgcaagggtggaactatcgggtatcacctgggtgttaactgcaagatcaagtcctgtactactcgtgcttctgtgactccaagaacgagtgctctggaccgacatgctctccaatttcggttacctggctaccagtcacaaac
tacgctgcatccggcagaaggcggtactgacgtggtaccgaggtggccccccggataaaagcatcatcaatgccacagacccctga