

PRKRIP1-PRKR interacting protein 1 (IL11 inducible) Gene

Reference: PTXBC014298

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C114; KRBOX3; PRKR-interacting protein 1; KRAB box domain containing 3; likely ortholog of mouse C114 dsRNA-binding protein; PRKR interacting protein 1 (IL11 inducible)		

For research use only.

GENE DESCRIPTION	PRKR interacting protein 1 (IL11 inducible)
GENE ID	79706
ORIGIN SPECIES	Human
ACCESSIONS	BC014298
VECTOR	pDONR223
NCBI SYMBOL	PRKRIP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctagccagccgctcctcgtgacaccagggcccaagaaagagccgagacgctcgtatcccaagaatcgccggaggagcagaagctcaagctggagcggctcatgaagaacccggacaaagcagttcattccagagaaatgagtgaatggcacctcgac
ctccccagaattgtccgagatgtcatgggttcaagtgctggggccggcagtgagagttccacgtgtacagacatctgcgccggagagaatcagcagacaggactacatggatgcatggctgagaagcaaaattgtatgcagagtttcagaaaagactgaaaagaataaattgc
tcagaggagcagaccgaaagcgggaagaagcgcagaagttaaaagagaagaattactggcaagaagatgaacacagaagaacaaaggacccggtcagcccaaggagcaggggtccagcagctctcggaggcatctggaacagaggaggaggaggaa
gtgccagtttcacatgggcatga