

CTNNBIP1-catenin, beta interacting protein 1 Gene

Reference: PTXBC014300

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	beta-catenin-interacting protein 1; beta-catenin-interacting protein ICAT; inhibitor of beta-catenin and Tcf-4; catenin beta interacting protein 1		

For research use only.

GENE DESCRIPTION	catenin, beta interacting protein 1
GENE ID	56998
ORIGIN SPECIES	Human
ACCESSIONS	BC014300
VECTOR	pDONR223
NCBI SYMBOL	CTNNBIP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaccgcgagggagctccgggaagagtcggaggagatgtacattcagcagaaggtccgagtgctgctcatgctcggaagatgggatcaaacctgacagccagcagaggaggagttctgcgcacctatgcaggggtgtcaacagccagctcagccagctgcctccgactc
catcgaccaggggtgcagaggacgtggtgatggcgtttccaggtcgagacggaagaccggaggcagtag