

ERH-enhancer of rudimentary homolog (Drosophila) Gene

Reference: PTXBC014301

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DROER; enhancer of rudimentary homolog; enhancer of rudimentary homolog (Drosophila)		

For research use only.

GENE DESCRIPTION	enhancer of rudimentary homolog (Drosophila)
GENE ID	2079
ORIGIN SPECIES	Human
ACCESSIONS	BC014301
VECTOR	pDONR223
NCBI SYMBOL	ERH
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctcaccatttgcgtgacagcctaccaagagccagaagcagaacttatctgactacgaatctgtgaatgaatgcagggtgtttgtaaatgtatgaagaacatctgaaaagaatgaatcccaacagtcctctatcacatgacatcagtcagttgtttgatttcacgatg
atctggcagaccctcagctgcctggtttaccgagctgataccagacataccagccttatacaaaagactggattaaagagaagatctacgtgctccttcgtcggcaggcccaacaggctgggaataa