

RPL10L-ribosomal protein L10-like Gene

Reference: PTXBC014310

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RPL10_5_1358; 60S ribosomal protein L10-like; ribosomal protein L10 like		
For research use only.			

GENE DESCRIPTION	ribosomal protein L10-like
GENE ID	140801
ORIGIN SPECIES	Human
ACCESSIONS	BC014310
VECTOR	pDONR223
NCBI SYMBOL	RPL10L
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggcgcgcgcagctcgtgttaccggtattgtaagaacaagccgtaccacaaatctgttctgccgaggggtcctgatccaagatccgcacatttgacctgggtagaagaaggcaaaagtggatgagttccactcgggtgccacatgggtctgatgaatatgagcagctgtctt
ctgaagccctggagggccgcccgtatttggccaacaatacatggfgaaggtgtggcagagatggctttcacatgcgagtgccgctccatccctccatgtcatccgcacacaagatgttctctgtctggggctgacaggtccagacaggtatgcgaggtgccttggaaaacccc
aggggtactgtagccgggtccacattggtcaagtcacatgtccatccgcaccaagcttcagaacgaggagcatgtgattgaagccttgcgcagggccaaagtcaagttccctggacgccagaagattcatatccaagaagtggggcttcacgaagttaatgctgacgaattgaagaca
tggtggccaagaagtgcctcatccctgatggtgtgagtcagtagtccagtcagtcgcccctggacaagtgccgggtctgcactcatga