

SNAPC5-small nuclear RNA activating complex, polypeptide 5, 19kDa Gene

Reference: PTXBC014315

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	snRNA-activating protein complex subunit 5; SNAPc 19 kDa subunit; SNAPc subunit 5; snRNA-activating protein complex 19 kDa subunit; small nuclear RNA activating complex polypeptide 5		

For research use only.

GENE DESCRIPTION	small nuclear RNA activating complex, polypeptide 5, 19kDa
GENE ID	10302
ORIGIN SPECIES	Human
ACCESSIONS	BC014315
VECTOR	pDONR223
NCBI SYMBOL	SNAPC5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctgagccggtcaggaactgcgaaggaggaggagacgtgctgcggtgaaggcagccctgcacgaccagctgaaccgcctcaagatgttggtcatgtagacaatgaagcatcaatcaaccaaacacccctggagctgagcacaagagtcgtgacggaagaggagga
ggaggaagaggagaagaatcagattcctaa