

PHLDA3-pleckstrin homology-like domain, family A, member 3 Gene

Reference: PTXBC014390

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TIH1; pleckstrin homology-like domain family A member 3; TDAG51/Ipl homolog 1; pleckstrin homology-like domain, family A, member 2; pleckstrin homology like domain family A member 3		

For research use only.

GENE DESCRIPTION	pleckstrin homology-like domain, family A, member 3
GENE ID	23612
ORIGIN SPECIES	Human
ACCESSIONS	BC014390
VECTOR	pDONR223
NCBI SYMBOL	PHLDA3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacggcgccggcgacggcctaccgtgctcaaggaggcgctgctggagaagcgacggcgccggcgctgctgagctgtggaaagcggaagcgtgctgctcaccgaacggcgccggcgctgagccaaggcgccggcgccggccaaggagctcagcttgcgccgca
tcaaggccgtggagtgctggagagcaccggcgccacatctacttcacgctggtgaccgaaggggcgccggcgagatcgacttccgctgccccctgggaagatccggctggaaacgccagatcacctaggctggtcaagttcaagaaccagcaggccatccagacagtgcgggc
ccggcagagcctcgggacgggacccctcgtgtcctaa