

MLANA-melan-A Gene

Reference: PTXBC014423

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MART-1; MART1; melanoma antigen recognized by T-cells 1; antigen LB39-AA; antigen SK29-AA; protein Melan-A; melan-A		

For research use only.

GENE DESCRIPTION	melan-A
GENE ID	2315
ORIGIN SPECIES	Human
ACCESSIONS	BC014423
VECTOR	pDONR223
NCBI SYMBOL	MLANA
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccaaagagaagatgctcacttcacatggttacccaagaaggggcacggccactcttacaccacggctgaagagccgctgggacgcacctgacagtgatcctgggagctctactgctcagcggctgttggtattgagaagacgaaatggatacagagccttgatgataaaagtcttcattgtggcactcaatgtgccttaacaagaagatgccacaagaaggggttgatcatcgggacagcaaaagtgtctcttcaagagaaaaactgtgaacctgtgttcccaatgctccacctgcttatgagaactctctgcagaacagtcaccaccactttacccttaa