

CD160-CD160 molecule Gene

Reference: PTXBC014465

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD160 molecule; CD160-delta Ig; CD160 transmembrane isoform; CD160 antigen; NK1; NK28; natural killer cell receptor BY55; natural killer cell receptor, immunoglobulin superfamily member		

For research use only.

GENE DESCRIPTION	CD160 molecule
GENE ID	11126
ORIGIN SPECIES	Human
ACCESSIONS	BC014465
VECTOR	pDONR223
NCBI SYMBOL	CD160
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctgttgaacccggcagaggtgctgtgccctggccatcctgctggcaattgtggacatccagtcgtggtgagtcattaacatcaccagtcagctccagggaaggacgcgactaaactaatctgtactgtatggcataagaagaaggctgaggggtttagtgtttttgtgca
aggacaggtctgagactgttccctgagaccagttaaacaagctgagacttaaaaggatcctgggatagatggtgttggatgaatatacatctcagttgattggtcaccataagccaagtcacaccgttgacagtgaggacctaccagtgtgtgccagaagccagaagtcaggtatccgcct
tcagggccatttttctccattctattcacagagacagggaactacacagtgacgggattgaaacaaagacaacaccttgagttcagccataatgaaggcactctcagttcaggcttcctacaagaaaggctcgggtatgctggtcaccagcctgtgtgcccttcaagctttaa