



Reference: PTXBC014507

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CDCD1; CDDC; CMD1A; CMT2B1; EMD2; FPL; FPLD; FPLD2; HGPS; IDC; LDP1; LFP; LGMD1B; LMN1; LMNC; LMNL1; MADA; PRO1; lamin; 70 kDa lamin; lamin A/C-like 1; mandibuloacral dysplasia type A; prelamin-A/C; renal carcinoma antigen NY-REN-32		

For research use only.

GENE DESCRIPTION	lamin A/C
GENE ID	4000
ORIGIN SPECIES	Human
ACCESSIONS	BC014507
VECTOR	pDONR223
NCBI SYMBOL	LMNA
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]